

**Shower seat**  
Catalogue number: VCTP0021

**VITEA®**  
**CARE**



**Product characteristics:**

- aluminium structure,
- seat and back made of easy cleanable polyurethane,
- equipped with handles,
- flow openings,
- step regulation of height every 2.5 cm.

**Use:**

Professional supporting equipment, simplifying / enabling the patient to shower.

Min. Height:	36,0 cm
Max. Height:	46,0 cm
Width:	48,0 cm
Depth:	44,0 cm
Weight:	2,5 kg
Maximum weight support:	90,0 kg