



Small shower stool
Catalogue number: VCTP0011

VITEA[®]
CARE

Product characteristics:

- aluminium structure,
- made of easy cleanable polyurethane,
- equipped with handles,
- flow openings,
- step regulation of height every 2.5 cm.

Use:

Professional supporting equipment, simplifying / enabling the patient to shower.

| | |
|----------------------|----------|
| Min. Height: | 36,0 cm |
| Max. Height: | 46,0 cm |
| Width: | 48,0 cm |
| Depth: | 44,0 cm |
| Weight: | 2,0 kg |
| Max. Weight support: | 100,0 kg |

