



## UMBRELLA

children buggy

Catalogue number: DRVG0C

**VITEA**<sup>®</sup>  
**CARE**

### Product characteristics:

- aluminium, foldable Frome
- polyester fabric seat
- polyester backrest of adjustable inclination (two positions)
- aluminium footrests of adjustable length and foot stabilizing belts
- swivel front wheels with blockade
- swivel rear wheels with shock absorber and brake
- front guard and seat belt to protect the child from falling out
- crotch strap
- standard equipment: canopy, sleeping bag, rain cover, headrest, shopping basket, bag.

### Use:

For children and teenagers with cerebral palsy, spina bifida and muscle disorders.

Footrest height: 20 – 35 cm

Seat width: 39 cm

Seat depth: 30 cm

Backrest height: 54 cm

Backrest tilt angle: 0 – 40 °

Max stroller width: 60 cm

Max stroller length: 100 cm

Max stroller height: 116 cm

Front wheels diameter: 18,5 cm

Rear wheels diameter: 22,5 cm

Wheelchair weight: 11,5 kg

Dimension while folder: 126X50X40 cm

Max child weight: 35 kg

Max child height: 130 cm