

Aluminium wheelchair Active Sport
Catalogue number: VCWK9AS



VITEA[®]
CARE

 42/46 cm	 14,5 kg
 43 cm	 max 115 kg
 dt. 100 sz. 71 w. 96	

Product characteristics:

- light, foldable, aluminium frame of trestle construction,
- black, detachable nylon upholstery,
- seat and backrest based on Velcro bands which allow depth adjustment,
- foldable backrest,
- corrugated anti-bedsore cushion,
- safety hip belt to prevent falling out,
- swing-away, detachable footrests equipped with heel strap,
- regulated footrest length,
- regulated height of wheels position towards the Frame,
- rear anti-flip wheels of regulated height,
- full front wheel – size: 6",
- pneumatic rear wheels with quick disassembly system, size: 24",

- inflator for inflating rear wheels,
- brakes with pressure adjustment,
- aluminium rims,
- modern construction and design,
- frame colour: silver, red, blue.

Use:

General use comfortable and versatile wheelchair for people with reduced mobility, to use at home and for strolls.