

Aluminium wheelchair Comfort

Catalogue number: VCWK9C

VITEA[®]
CARE

 20,5 kg	 46 cm
 max 115 kg	 43 cm
 dt. 109,0 sz. 65,0 w. 103,5	



Product characteristics:

- foldable, aluminium frame of trestle construction,
- detachable nylon black upholstery,
- backrest based on Velcro bands which allow depth adjustment,
- seat made of nylon, equipped with thick, soft cushion,
- safety hip belt to prevent falling out,
- swing-away, detachable footrests with tilt angle regulation and calf pads,
- regulated footrest length,
- possibility of moving the centre of gravity due to adaptable position of rear wheel axle,
- swing-away detachable armrests,
- regulated handle height,
- rear barrier which provides additional stabilization and facilitates folding of the wheelchair,
- rear anti-flip wheels of regulated height,
- full front wheels, size 200/50,
- pneumatic front wheels equipped with quick dis-assembly system, size 24",
- inflator for inflating rear wheels,
- cane holder,
- brakes with pressure adjustment,
- aluminium rims,
- modern construction and design,
- frame colour: silver.

Use:

General use comfortable and versatile wheelchair for people with reduced mobility, to use at home and for strolls.