

**Special wheelchair stabilising back
and head Comfort Plus**
Catalogue number: VCWK7CP

VITEA[®]
CARE



Product characteristics:

- stable, reinforced steel frame,
- foldable backrest,
- detachable, anti-bedsore upholstery made of 'memory foam', equipped with waterproof cover,
- seat equipped with soft, anti-bedsore cushion made of 'memory foam' covered with waterproof and easy to clean fabric,
- anatomically molded backrest and seat,
- smooth regulation of backrest height,
- smooth regulation of seat tilt angle (together with backrest),
- smooth regulation of backrest tilt angle,
- regulated footrest height,
- swing-away detachable footrests equipped with smooth tilt angle regulation and calf pads,
- regulated height of front wheels axle,
- detachable armrests of adjustable height made of soft, comfortable and easy to wash fabric,
- smooth regulation of armrest horizontal position,
- regulated seat width,
- double handle of adjustable height,
- lateral wedges of adjustable horizontal and vertical position,
- rear anti-flip wheels of regulated height,
- detachable headrest of regulated height and tilt angle,
- brakes with pressure adjustment,
- aluminium rims,
- full front wheels, size 200/50,
- pneumatic rear wheels equipped with quick disassembly system, size 24",
- inflator for inflating rear wheels,
- modern and interesting design,
- frame colour: silver.

Use:

General use comfortable and versatile wheelchair for people with reduced mobility, to use at home and for strolls.