

## Special wheelchair stabilising back and head with toilet function

Catalogue number: VCWK7T

**VITEA<sup>®</sup>**  
**CARE**



### Product characteristics

- recliner type wheelchair, with the possibility of sloping the back from sitting to lying position,
- foldable aluminium frame,
- black seat made of polyurethane material,
- sanitary container,
- small table,
- removable and lifted footrests with padded support for calves, and regulated tilt angle,
- removable armrests,
- anti-overtum system,
- front wheels full or pneumatic,
- rear wheels pneumatic or spoke,
- wheel lock.

### Use:

General use wheelchair, to use at home and for strolls with toilet function.