



Steel and aluminium wheelchair
Catalogue number: VCWK9AC

VITEA[®]
CARE

 42/46 cm	 50 cm	 max 100 kg
 44 cm	 25 cm	
 di. 108 ax. 62 w. 92	 13 kg	
	 32 cm	

Product characteristics:

- foldable aluminium frame,
- removable black and red, or black and blue upholstery made of sweat absorbing material,
- belts made of Velcro[®] band allowing for the back depth to be regulated,
- anti-bedsore cushion,
- safety belts for the user,
- double brake system (for patient and their companion),
- back folding mechanism situated next to armrests, allowing the reduction in the size of the wheelchair for transport,
- lifted armrests,
- lifted and removable footrests with calves securing belt,

- pneumatic rear wheels with a quick disassembly system,
- front wheels pneumatic or PU,
- cross structure – foldable (double cross),
- aluminium pulls,
- frame colours: platinum, blue and black.

Use:

General use wheelchair, at home and for strolls.