

**Steel wheelchair**  
Catalogue number: **VCWK4**

**VITEA<sup>®</sup>**  
**CARE**

 42/45 cm	 44,5 cm
 31 cm	 H. 107 sz. 62/65 w. 89
 18 kg	 max 100 kg



**Product characteristics:**

- foldable, steel frame,
- nylon black upholstery,
- anti-bedsore cushion,
- lifted armrests,
- lifted and removable footrests,
- belt to prevent falling out,
- pneumatic rear wheel 24" (additionally equipped with quick disassembly system),
- full front wheels 200X50,
- cross structure – foldable (double cross),
- rear pocket,
- frame colour: platinum.

**Use:**

General use wheelchair, at home and for strolls.